

Getting Ready to Restart Your Lab? We Can Help

Recommended GC and GC/MS consumables for bringing your lab back online

Are you looking to open your doors, turn on your instruments, and restart your lab? Agilent is here to make sure you have everything you need, every step of the way.

With that in mind, Agilent technical experts have compiled a list of recommended products designed to:

- Quickly get your operations running as efficiently and productively as possible—particularly if your lab has been shut down for multiple weeks.
- Minimize the risk to your GC flow path.
- Simplify installation and get your lab back to producing reliable results.

Designed for exceptional quality, performance, and lifetime, these products allow you to spend more time analyzing samples, and less time worrying about instrument downtime.

It's all part of our commitment to helping your lab navigate its "new normal" after these unprecedented times.



Easy ordering information

Want to add the items under each category to your My Favorites list at the Agilent online store? Simply click the header links. Then, enter the quantities for the products you need. Your list of items will remain under My Favorites for you to use with future orders.

Airbath GC Consumables	Part Number
Collared column nut, self-tightening inlet/detector	G3440-81011
Replacement collar, for self-tightening nut	G3440-81012
Splitless, UI, low-frit liner, 4 mm	5190-5112
Splitless, UI, low-frit liner, 4 mm, 5/pk	5190-5112-005
Universal, UI, mid-frit liner, 4 mm	5190-5105
Universal, UI, mid-frit liner, 4 mm, 5/pk	5190-5105-005
Cap insert, for diffusion caps, for 4 mL vials, 12/pk	07673-40180
FID jet, universal fit, 0.011 inch id	5200-0176
FID jet, universal fit, 0.018 inch id	5200-0177
Split vent trap PM kit, single cartridge	5188-6495
Septa, nonstick BTO inlet, 11 mm, 50/pk	5183-4757
Septa, nonstick adv green, 11 mm, 50/pk	5183-4759
Syringe, 10 µL tapered, FN 23-26s/42/HP	5181-1267
Ferrule, graphite, short, 0.5 mm id, 10/pk	5080-8853
Ferrule, 0.4 mm id, 15% graphite/85% Vespel, 10/pk	5181-3323
GC inlet seal, gold plated, with washer, 10/pk	5190-2209

Intuvo Consumables	Part Number
S/SL Guard Chip, 2/pk	G4587-60565
MMI Guard Chip, 2/pk	G4587-60665
S/SL Jumper Chip	G4587-60575
MMI Jumper Chip	G4587-60675
S/SL compression bolt	G4582-20085
Guard Chip compression bolt	G4581-20006
Compression bolt (quantity of 2)	G4581-60260
Polyimide gasket, 5/pk	5190-9072
Nickel gasket, 2/pk	5190-9073
Inlet chip	G4581-60031
Swaged MS tail	G4590-60009

GC/MS Consumables and Accessories	Part Number
Collared column nut, self-tightening MSD	G3440-81013
IDP-3 scroll pump for Turbo 597X	G6696A
IDP-10 scroll pump for GC/TQ	G6697A

Sample Containment	Part Number
Vial, screw top, clear, write-on spot, certified, 2 mL, 100/pk	5182-0715
Vial, screw top, amber, write-on spot, certified, 2 mL, 100/pk	5182-0716
Vial, crimp top, clear, write-on spot, certified, 2 mL, 100/pk	5182-0543
Vial, crimp top, amber, write-on spot, certified, 2 mL, 100/pk	5181-3376
Cap, screw, blue, PTFE/red silicone septa, 100/pk	5182-0717
Cap, crimp, silver aluminum, PTFE/silicone septa, 100/pk	5182-0552

GC Handheld Diagnostic Tools	Part Number
ADM Flow Meter	G6691A
Replacement cartridge for ADM Flow Meter	G6692A
Electronic leak detector	G3388B

Replacement Gas Purifiers	Part Number
Gas Clean filter, oxygen	CP17970
Gas Clean filter, moisture	CP17971
Gas Clean filter, hydrocarbon	CP17972
Gas Clean filter, carrier gas	CP17973
Big universal trap, 1/8-inch fittings, helium	RMSH-2
Big universal trap, 1/4-inch fittings, helium	RMSH-4
Big universal trap, 1/8-inch fittings, hydrogen	RMSHY-2
Big universal trap, 1/4-inch fittings, hydrogen	RMSHY-4
Big universal trap, 1/8-inch fittings, nitrogen	RMSN-2
Big universal trap, 1/4-inch fittings, nitrogen	RMSN-4

www.agilent.com/chem/gc-supplies

DE.524039519

This information is subject to change without notice.

© Agilent Technologies, Inc. 2020
Printed in the USA, June 10, 2020
5994-2099EN